

Overview WiKo anterior clavicle plates 3.5

titanium, TiOB[®]-coated



pict.: 5.3802.06



pict.: 5.3812.08

Item no.	Item description
----------	------------------

Anterior lateral clavicle plates

5.3812.08	WiKo anterior lateral clavicle plate 3.5; 8 holes; length 105 mm
5.3812.09	WiKo anterior lateral clavicle plate 3.5; 9 holes; length 117 mm
5.3812.10	WiKo anterior lateral clavicle plate 3.5; 10 holes; length 128 mm

Anterior medial clavicle plates

5.3802.06	WiKo anterior medial clavicle plate 3.5; 6 holes; length 80 mm
5.3802.07	WiKo anterior medial clavicle plate 3.5; 7 holes; length 91 mm
5.3802.08	WiKo anterior medial clavicle plate 3.5; 8 holes; length 103 mm

Information

In Germany the system name WiKo[®] is a Registered Trade Mark. The system name WiKo[®] is derived from the German designation "Winkelstabile Kombination", which translates into English as "Angle-stable Combination". This is a plating technology which offers stable locking at a specified angle. It includes all the angle-stable plates both with and without compression.

In Germany, in the EU states, Switzerland and in the USA TiOB[®] is a Registered Trade Mark.

We hope, this offer meets your expectations.
We would be happy to receive your order.



28/01/2015/KaMo

Königsee Implantate GmbH
OT Aschau | Am Sand 4
07426 Allendorf | Germany
Fon + 49 36738 498-360
Fax + 49 36738 498-559
E-Mail: export@koenigsee-implantate.de

For further information please visit:
www.koenigsee-implantate.com