

Set angle-stable distal and proximal Radius plates

for screws D 2,0/3,5 mm



Picture: Radial plate palmar, right, titanium

Code N°: 5.898.63

In head: 6 holes for the angle-stable screw connection (screws D 2,0 mm with conic thread in head are placed in the picture)

In shank: 3 shank holes for cortical screws D 3,5 mm




Picture: Radial plate dorsal, right, titanium

Code N°: 5.892.53

In head: 5 holes for the angle-stable screw connection (screws D 2,0 mm with conic thread in head are placed in the picture)

In shank: 3 shank holes for cortical screws D 3,5 mm



Code N°
Set complete
Titanium
19.436.00



Picture: Proximal Radius plate for screws D 2,0 mm, Titanium
Code N°: 5.884.23

■ Angle-stable distal dorsal Radius plate, Titanium

Description	Code N°	Quantity
Right, for screws D 2,0/3,5 mm, 63x10 mm	5.892.53	1
Left, for screws D 2,0/3,5 mm, 63x10 mm	5.893.53	1

■ Angle-stable distal palmar Radius plate, Titanium

Description	Code N°	Quantity
Right, for screws D 2,0/3,5 mm, 50x10 mm	5.898.63	1
Left, for screws D 2,0/3,5 mm, 50x10 mm	5.899.63	1

■ Angle-stable proximal Radius plate, Titanium

Description	Code N°	Quantity
2/3 holes, for screws D 2,0 mm, 32x7 mm	5.884.23	1
2/4 holes, for screws D 2,0 mm, 40x7 mm	5.884.24	1

■ Instruments

Description	Code N°	Quantity
Drill bit for quick coupling D 1,5 mm / L 85 mm	2.904.20	1
Drill bit for quick coupling D 2,5 mm / L 105 mm	2.904.05	1
Drill guide for angle stable screws 2.0 mm	2.977.06	2
Drill guide for angle stable screws 3,5 mm	10.500.27	1
Hex screwdriver with handle, wrench width 1.5 mm	2.940.15A	1
Hex screwdriver with handle, wrench width 2.5 mm	2.940.25SA	1
Screw gauge for screws with conical head thread	2.953.51	1
Screw forceps, self holding	2.954.01	1
Perforated autoclavable container	19.435.00	1

■ Wire

Description	Code N°	Anzahl
Kirschner wire with trocar point/ round end, Ø1.4x150 mm	6.031.14	5

■ Angle-stable cortical screw, Titanium

D 2,0 mm (fully threaded, self-tapping)

Length	Code N°	Quantity
16 mm	3.114.16	5
18 mm	3.114.18	5
20 mm	3.114.20	5
22 mm	3.114.22	5
24 mm	3.114.24	5
26 mm	3.114.26	5
28 mm	3.114.28	5
30 mm	3.114.30	5



■ Cortical screw, Titanium

D 2,0 mm (fully threaded, self-tapping)

Length	Code N°	Quantity
14 mm	3.110.14	5
16 mm	3.110.16	5
18 mm	3.110.18	5
20 mm	3.110.20	5
22 mm	3.110.22	5
24 mm	3.110.24	5



■ Cortical screw, Titanium

D 3,5 mm (fully threaded, self-tapping)

Länge	Bestell Nr.	Anzahl
14 mm	3.132.14	5
16 mm	3.132.16	5
18 mm	3.132.18	5
20 mm	3.132.20	5
22 mm	3.132.22	5
24 mm	3.132.24	5



KÖNIGSEE IMPLANTATE UND INSTRUMENTE

ZUR OSTEOSYNTHESE GmbH

AM SAND

D-07426 ALLENDORF

Phone +49 (0) 36738 498-0

Fax +49 (0) 36738 498-19

e-Mail: info@koenigsee-implantate.de

Set angle-stable distal and proximal Radius plates

for screws D 2,0/3,5 mm

